



**NewBorn Paste**  
from  
**KeyAg**  
DISTRIBUTORS

## Excell PRO NewBorn Paste

KeyAg utilizes Egg Protein Technology (EPT) in our **Excell PRO** line up of products. EPT can produce specific IgY immunoglobulins in concentrations that are clinically proven to be effective in neutralizing intestinal pathogens. EPT specific egg proteins can also aid in the passive transfer of critical immune factors (IgY).

These concentrations are determined through a patented process called the ELISA Technique. ELISA (Enzyme Linked Immunosorbent Assay) is highly accurate and the most economical way to determine these levels.

Determining the concentration of specific egg proteins that are produced through EPT is critical to formulating effective end products like the **Excell PRO** family of products.

### With six specific IgY Immunoglobulins



**Use Birth to 3 Months**

**Contains Procrypt**

**Proudly manufactured in the USA by:**

**KeyAg**  
DISTRIBUTORS  
[www.keyag.com](http://www.keyag.com)



**4525 E, 3425 N**  
**Murtaugh, Idaho 83344**  
**1.800.388.3659**  
**info@keyag.com**